



DATA SHEET

EpisoZYme Taq DNA Polymerase with Standard Taq Buffer

Cat. No.	Pack Size	Conc.
EZ001S	500U	5U/ μ L
EZ001M	1000U	5U/ μ L
EZ001L	5000U	5U/ μ L

EpisoZYme Taq DNA Polymerase is a recombinant DNA polymerase purified from *Escherichi coli* that carries the Taq DNA Polymerase gene of *Thermus aquaticus*. It has 5'-3' DNA polymerase activity and 5'-3' exonuclease activity.

Reagents Provided

- **Taq DNA Polymerase in Storage Buffer** (100 mM KCl, 10 mM Tris-HCl: pH 7.4, 0.1 mM EDTA, 1 mM DTT, 0.5% Tween 20, and 50% glycerol)
- **10x Reaction Buffer:** 100 mM Tris-HCl (pH 8.8 at 25°C), 500 mM KCl, 1% Triton X-100
- **25 mM MgCl₂ Solution**

Recommended PCR Reaction Mix

10x Buffer	2.5 μ l	5 μ l	1x
2.5 mM dNTPs	2 μ l	4 μ l	200 μ M
10 μ M For. Primer	1.25 μ l	2.5 μ l	0.5 μ M (0.05–1 μ M)
10 μ M Rev. Primer	1.25 μ l	2.5 μ l	0.5 μ M (0.05–1 μ M)
Taq DNA Poly. (5U/ μ l)	0.1 μ l	0.2 μ l	1 units / 50 μ l PCR
25 mM MgCl ₂	2 μ l	4 μ l	2 mM (1.5-2.5 mM)
Template DNA	variable	variable	<100 ng
Nuclease-Free Water	to 25 μ l	to 50 μ l	

Application


- Primer Extension
- Colony PCR
- Microarray Analysis

Thermocycling Conditions for a Routine PCR

Step	Temp.	Time
Initial Denaturation	95 °C	1 m
25-30 Cycles	95 °C 45-60 °C 74 °C	15 s 10-30 s 20 s/kb
Final Extension	74 °C	5 m
Hold	r. t.	

Quality Control Assays

PCR; cDNA cloned into a plasmid, 30 cycles of PCR amplification of 5 ng template DNA with Taq DNA Polymerase in the presence of 200 μ M dNTPs and 0.2 μ M primers in standard Taq reaction buffer results in the expected 3 kbp product.



Safety warnings and precautions

This product is designed for research purposes and *in vitro* use only. According to common laboratory safety practice, it is recommended to wear protective clothing, gloves and safety glasses.

